

# Swiftpet Pro

## Pipet Controller

Pipetting using glass and plastic volumetric pipets has never been easier than with the Swiftpet Pro Pipet Controller. An ergonomically shaped handle along with optimally located, smooth pushbuttons and switches guarantee effortless pipetting even during periods of extended use. A convenient charging stand allows the pipet controller to be stored when not in use. Easily accessible switches on the Swiftpet Pro allow the user to select different operation modes depending on the volume of pipet and viscosity of liquid. The Swiftpet Pro is an excellent tool for pipetting in tissue culture laboratories, water treatment facilities, academic labs, and more.

### Safe and efficient work

- ▶ Protected by a PTFE filter that blocks liquids from entering the unit
- ▶ Autoclavable filter, pipet holder and nosepiece
- ▶ UV-resistant body for safe sterilization
- ▶ Powerful, environmentally friendly 3 Ni-MH batteries enable many hours of continuous work



### Comfort while working

- ▶ Suitable for glass and plastic volumetric pipets 0.5 - 100 mL
- ▶ Ergonomically shaped handle
- ▶ Optimally located function buttons
- ▶ Convenient charging stand for bench/wall mounting
- ▶ Easy access to batteries



- ▶ LCD display showing battery charge level



### Adjustment of speed and working mode

- ▶ Easily accessible function buttons for control of speed and working mode
- ▶ Additional speed adjustment via the pressure applied to the trigger buttons



Product	Cat. No.
Swiftpet Pro with charging stand	0390
Charging stand	SP29056
Pipettor holder	SP29054
Nosepiece	SP29053
Battery cover with "wings"	SP29064
Micro filter 0.20 µm, pack of 5	SP9143
Micro filter 0.45 µm, pack of 5	SP9144
Universal charger 9V type: EU, US, UK, AU	SP29100
Hanger	SP9029