

| 型 號     | FDM-2                        | FDM-5      | FDM-12     | FDM-20     |
|---------|------------------------------|------------|------------|------------|
| 冷凝器溫度   | 啟動後10分鐘內可降至-50/-84℃以下(室溫20℃) |            |            |            |
| 最低真空度   | 幫浦啟動後10分鐘內可降至40mTorr以下       |            |            |            |
| 水凍結量    | 2 L                          | 5 L        | 12 L       | 20 L       |
| 歧管數     | 8                            | 12         | 16         | 24         |
| 熱氣融冰    | Option                       | Yes        | Yes        | Yes        |
| 電壓      | 110/220V                     | 220V       | 220V       | 220V       |
| 尺寸W×D×H | 46×50×56cm                   | 50×60×78cm | 82×75×80cm | 82×75×80cm |



# 歧管式冷凍乾燥機

## Manifold Freeze Dryer



| 冷凍乾燥瓶      |                |              |           |         |
|------------|----------------|--------------|-----------|---------|
| Flask Size | Complete Flask | Flask Bottom | Flask Top | Filter  |
| 150ml      | 1010150        | 1000150      | 1012000   | 1012001 |
| 300ml      | 1010300        | 1000300      | 1012000   | 1012001 |
| 600ml      | 1010600        | 1000600      | 1011000   | 1011001 |
| 1200ml     | 1011200        | 1001200      | 1011000   | 1011001 |

| 真空幫浦           |                                                                                  |  |          |
|----------------|----------------------------------------------------------------------------------|--|----------|
| Catalog Number | Description                                                                      |  | Oil Mist |
| 1050205        | RV5 , Pumping speed 103 L/min , Ultimate total pressure2*10 <sup>-3</sup> mbar   |  | 1065005  |
| 1050136        | GLD-N136 , Pumping speed 165 L/min , Ultimate pressure 6.7*10 <sup>-3</sup> mbar |  | 1065012  |
| 1050212        | RV12 , Pumping speed 283 L/min , Ultimate total pressure2*10 <sup>-3</sup> mbar  |  | 1065012  |

| Catalog Number | Description    |
|----------------|----------------|
| 1051001        | 1L高真空幫浦油       |
| 1051004        | 4L高真空幫浦油       |
| 1080001        | 手動封瓶裝置         |
| 1030001        | 推車(FDM-2)      |
|                | IQ/OQ Document |



智慧型全自動控制  
觸控式彩色螢幕

